

## Product Information Sheet

<b>Catalog Number:</b>	A1009
<b>Product Name:</b>	Anti-huB7-H3 mAb
<b>Product Description:</b>	Recombinant human IgG1 antibody against human B7-H3
<b>Clone Name:</b>	A84
<b>Protein Concentration:</b>	1.0 mg/mL
<b>Protein Formulation:</b>	Sterile 1X PBS, pH 6.8
<b>Storage Conditions:</b>	-20°C or below, avoid repeated freeze and thaw
<b>Protein Stability:</b>	>30 days if stored at 4°C or below
<b>Estimated Purity:</b>	~95%
<b>Protein Aggregation:</b>	Not obvious by SDS-PAGE
<b>Specificity Analysis:</b>	Binding of anti-human B7-H3 mAb A84 (Cat.# A1009) to CHO cells expressing full length human B7-H3.

