

Product Information Sheet

Catalog Number:	A1012
Product Name:	Anti-human CD47 huIgG4
Product Description:	Recombinant human IgG4 antibody against human CD47
Clone Name:	AB6
Protein Concentration:	1.0 mg/mL
Protein Formulation:	Sterile 1X PBS, pH 6.8
Storage Conditions:	-20 °C or below, avoid repeated freeze and thaw
Protein Stability:	>30 days if stored at 4°C or below
Estimated Purity:	~95%
Protein Aggregation:	Not obvious by SDS-PAGE
Specificity Analysis:	Anti-human CD47 huIgG4 mAb, clone AB6, binding to human CD47 expressed on CHO cells

